

# The New NanoScope V Controller:

## *New power and capabilities for MultiMode V, Dimension V, NanoMan V, and PicoForce*

### Fast, Dependable Data Capture

The new NanoScope® V controller utilizes advanced electronics, including A/D and D/A converters operating at 50MHz, to deliver reliable, high-speed data capture. This state-of-the-art controller allows researchers to record and analyze tip-sample interactions (e.g., at pull-off in force spectroscopy) that probe nanoscale events at time-scales previously inaccessible to SPM.

High speed data capture is simultaneous with imaging or ramping and is independent of microscope mode.

### Flexible Controller Features

The NanoScope V enables up to eight images to be simultaneously displayed in real-time (and captured for analysis)



The NanoScope V controller.

with unprecedented signal-to-noise ratio. The controller incorporates three independent lock-in amplifiers and provides thermal tune measurements of cantilever resonances up to 2MHz.

The controller also affords easy access to most input and output signals through front-panel BNCs. Input data into the controller from an external source (e.g., photomultiplier tube) is supported, as is user access to lock-in

amplifiers and to signals to/from a microscope (e.g., XYZ sensors, amplitude, phase).

### Highest Pixel Density

Up to 5000 x 5000 pixel images, improve “time to results” by reducing the need to capture several images at lower pixel densities; and eliminating the requirement for offset adjustments to correlate information from multiple images. The high pixel density also allows observation of large structures and small features in the same image.

### Outstanding Software Functionality

The new controller’s NanoScript open-architecture option provides an expansive set of SPM control functions for custom experiments and nanoscale research (e.g., nanomanipulation in X, Y, Z; automated scanning; nanolithography with different tip-sample interactions).

These functions can also be called from any programming client of Microsoft’s Component Object Model (COM), including LabView, MATLAB, Visual Basic, Ruby, Python, C++/MFC, Excel®, and Word®.

### EASY-AFM, Remarkable Simplicity

For the ultimate in streamlined operational simplicity, the NanoScope V controller’s Easy-AFM™, ease-of-use feature, offers an intuitive, easy-to-follow graphic user interface. Easy-AFM reduces the time for initial setup (including probe, laser, and detector alignments), engaging the sample with the probe, adjusting the scanning parameters, and obtaining TappingMode™ images on most samples.

### Highlights of Features and Benefits

- High-speed electronics enable recording and analyzing data (e.g., force pulling experiments) from nano-scale events at timescales previously inaccessible to SPM
- Advanced software features simplify operation for novice or infrequent users, and give experienced users more power, and more ways to extract information from data
- High pixel density images up to 5k x 5k massively improve productivity
- Up to 8 simultaneous channels in real-time-scanning and offline-visualization/analysis enable correlating information about unprecedented number of sample properties
- Open architecture provides new options to design and run customized experiments, including with third-party software
- Differential signals and adjustable gains improve signal-to-noise ratio
- Closed-loop XYZ scanner compatibility utilizes the latest in scanner design and functionality
- Expandable modular design is forward-looking towards supporting additional functions in future
- Customizable signal access using software-controlled signal routing expands options for experimentation

## Hardware

### Electronics

#### High speed and expanded AC capabilities

- Two high-speed (50MHz) 14-bit ADCs sample and digitize the probe signal (+/-2V)
- Two high-speed (50MHz) 16-bit DACs provide sinusoids (+/-10V), plus 16-bit offset bias superposed (+/-10V, adjustable up to 500kHz)
- Nine mid-speed (500kHz) 18-bit ADCs (+/-12V) enable multiple digital feedback loops to operate at 2 $\mu$ s speed, allowing faster scanning and data capture with fast actuators.
- Two independent high-speed lock-ins (1KHz – 5MHz) and one mid speed lock-in (0.1Hz – 50KHz) allow sampling and determining amplitude and phase of up to three independent signals, or analyzing higher harmonics of a signal (0.5KHz – 2.1MHz or 5Hz – 21KHz lock-in bandwidths)
  - Three independent reference signals available (2 at up to 5MHz, 1 at up to 50KHz)

#### User-accessible hardware input/output

- Three mid-speed ADCs, two high-speed ADCs
- Two digital inputs, two digital outputs (TTL compatible)
- Two high-speed and three mid-speed 16-bit DACs (+/-10V)

### Hardware options

- “Pico Angler” hand-held module for tactile interface and feedback in force spectroscopy
- Signal Access Module for “MultiMode<sup>®</sup> V”

## Software

- Exceptional Force Spectroscopy control and analysis
  - “Point and Shoot”: precisely located force or tunneling spectroscopy at any XY coordinate in an image, with a single click of the mouse
  - Force Volume imaging up to 256 x 256 x 256 data-points per image
- Open architecture using easy to use Applications Programming Interface (API) and Component Object Model (COM)
- High-speed data capture commensurate with high-speed electronics
- Easy-AFM, the “Expert System” graphic user interface, easy-to-follow for new or infrequent SPM users to set up (including probe, laser, and detector alignments); engage the sample with the probe in air; adjust the scanning parameters; and obtain high- quality, useful TappingMode images on most samples with minimal user intervention

- Expanded, enhanced image capture, visualization, and analysis
  - Flexible combinations of sample-per-line and lines-per-frame imaging (e.g., 512 x 512, 8192 x 3200, 16384 x 1600, 5000 x 5000)
  - Up to 8 channels of data can be simultaneously collected
  - Two- and three-dimensional image display options with 24 color tables
  - Full suite of image analysis at high pixel densities
  - Continuous data capture from successive scans of the same area
  - Automatic data capture for calibration
- Integrated graphic user interface for
  - Thermal tuning of AFM cantilevers with fundamental resonance frequencies up to 2MHz (e.g., for spring constant measurements)
  - Controlling lock-in hardware
  - Controlling cantilever Q
  - Adjusting scan rate, size, and offset (calibrations remain unchanged)
  - Disabling scanning along the slow scan axis
  - Adjusting gain, set point, drive frequency/amplitude/phase
  - Selecting AFM and STM feedback modes
  - Applying bias voltages to tip and/or sample
  - Adjusting Z drive voltage range (11V to 440V) commensurate with desired vertical resolution
  - Selecting gains for most input signals
  - Selecting and setting parameters for multiple auxiliary input channels
  - Correcting for XY linearity
- Optional: Nanomanipulation, NanoScript, TRmode

### Software Utilities

- ASCII, jpeg, BMP export
- Printing images, spectroscopy curves, cantilever tune plots, image analysis windows

## Other Specifications

|                                    |   |
|------------------------------------|---|
| <b>Computer</b>                    | Premium computer system<br><i>(call for current specifications)</i> |
| <b>Space/footprint requirement</b> | 20" x 30" footprint,<br>10" height clearance                        |
| <b>Power requirement</b>           | 520W peak power consumption   |
| <b>Weight</b>                      | 46 pounds   |

*Note: Performance specifications are typical and subject to change without notice.*



Veeco Instruments Inc.

112 Robin Hill Road  
Santa Barbara, CA 93117

805-967-1400 • 1-888-24-VEECO

Find out more at [www.veeco.com](http://www.veeco.com)

DS73, Rev A0

© 2006 Veeco Instruments Inc. All rights reserved.  
Easy-AFM is a trademark of Veeco Instruments Inc.  
NanoScope and MultiMode are registered trademarks of Veeco Instruments Inc.  
Excel and Word are registered trademarks of MicroSoft Corp.  
MATLAB is a registered trademark of The Mathworks, Inc.